

兔抗 LGR5 多克隆抗体

中文名称：兔抗 LGR5 多克隆抗体

英文名称：Anti-LGR5 rabbit polyclonal antibody

别名：FEX, HG38, GPR49, GPR67, GRP49

抗原：LGR5

储存：冷冻（-20℃）

宿主：Rabbit

相关类别：一抗

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes Orphan receptor. Stem cell marker of the intestinal epithelium and the hair follicle. Target gene of Wnt signaling.
Applications:	ELISA, IHC
Name of antibody:	LGR5
Immunogen:	Fusion protein of human LGR5
Full name:	leucine-rich repeat containing G protein-coupled receptor 5
Synonyms :	FEX, HG38, GPR49, GPR67, GRP49
SwissProt:	O75473
ELISA Recommended dilution:	1000-5000

IHC positive control:	Human brain
IHC Recommend dilution:	25-100

