

## 小鼠抗 ADK 单克隆抗体

中文名称： 小鼠抗 ADK 单克隆抗体

英文名称： Anti-ADK mouse monoclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃）

宿 主： Mouse

抗 原： ADK

反应种属： Human

标 记 物： Unconjugate

克隆类型： mouse monoclonal

### 技术规格

<b>Background:</b>	This gene an enzyme which catalyzes the transfer of the gamma-phosphate from ATP to adenosine, thereby serving as a regulator of concentrations of both extracellular adenosine and intracellular adenine nucleotides. Adenosine has widespread effects on the cardiovascular, nervous, respiratory, and immune systems and inhibitors of the enzyme could play an important pharmacological role in increasing intravascular adenosine concentrations and acting as anti-inflammatory agents. Multiple transcript variants encoding different isoforms have been found for this gene.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	ADK

<b>Immunogen:</b>	Fusion protein of human ADK
<b>Full name:</b>	adenosine kinase
<b>Synonyms:</b>	AK
<b>SwissProt:</b>	P55263
<b>WB Predicted band size:</b>	41 kDa
<b>WB Positive control:</b>	HEK293 whole cell lysate
<b>WB Recommended dilution:</b>	1000-5000

