

兔抗 LHCGR 多克隆抗体

中文名称：兔抗 LHCGR 多克隆抗体

英文名称：Anti-LHCGR rabbit polyclonal antibody

别名：HHG; LHR; LCGR; LGR2; ULG5; LHRHR; LSH-R; LH/CGR; LH/CG-R

相关类别：一抗

抗原：LHCGR

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes the receptor for both luteinizing hormone and choriogonadotropin. This receptor belongs to the G-protein coupled receptor 1 family, and its activity is mediated by G proteins which activate adenylate cyclase. Mutations in this gene result in disorders of male secondary sexual character development, including familial male precocious puberty, also known as testotoxicosis, hypogonadotropic hypogonadism, Leydig cell adenoma with precocious puberty, and male pseudohermaphroditism with Leydig cell hypoplasia.
Applications:	ELISA, IHC
Name of antibody:	LHCGR

Immunogen:	Synthetic peptide of human LHCGR
Full name:	luteinizing hormone/choriogonadotropin receptor
Synonyms:	HHG; LHR; LCGR; LGR2; ULG5; LHRHR; LSH-R; LH/CGR; LH/CG-R
SwissProt:	P22888
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer and human lung cancer
IHC Recommend dilution:	25-100



