

兔抗 LIG1 多克隆抗体

中文名称:	兔抗 LIG1 多克隆抗体
英文名称:	Anti-LIG1 rabbit polyclonal antibody
相关类别:	一抗
储 存:	冷冻 (-20℃)
抗 原:	LIG1
宿 主:	Rabbit
反应种属:	Human, Mouse, Rat
标记物:	Unconjugate
克隆类型:	rabbit polyclonal

技术规格

Background:	This gene encodes a member of the ATP-dependent DNA ligase protein family. The encoded protein functions in DNA replication, recombination, and the base excision repair process. Mutations in this gene that lead to DNA ligase I deficiency result in immunodeficiency and increased sensitivity to DNA-damaging agents. Disruption of this gene may also be associated with a variety of cancers. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	LIG1
Immunogen:	Fusion protein of human LIG1
Full name:	ligase I, DNA, ATP-dependent
SwissProt:	P18858

ELISA Recommended dilution:	2000-10000
IHC positive control:	Human esophagus cancer and Human breast cancer
IHC Recommend dilution:	100-300

