

兔抗 SPATA6 多克隆抗体

中文名称：兔抗 SPATA6 多克隆抗体

英文名称： Anti-SPATA6 rabbit polyclonal antibody

别名： HASH; SRF1; SRF-1

抗原： SPATA6

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The Spata6 gene consists of 15 exons ranging in size between 40 and 596 bp. It may play a role in spermatid maturation or sperm function. The Spata6 gene was localized to chromosome 5, region q34-35 in the rat and to chromosome 1, region p32-35 in the human.
Applications:	ELISA, IHC
Name of antibody:	SPATA6
Immunogen:	Fusion protein of human SPATA6
Full name:	spermatogenesis associated 6
Synonyms :	HASH; SRF1; SRF-1

SwissProt:	Q9NWH7
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer and human lung cancer
IHC Recommend dilution:	25-100

