

兔抗 SPATA7 多克隆抗体

中文名称：兔抗 SPATA7 多克隆抗体

英文名称： Anti-SPATA7 rabbit polyclonal antibody

别名： HSD3; LCA3; HSD-3.1; HEL-S-296

抗原： SPATA7

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene, originally isolated from testis, is also expressed in retina. Mutations in this gene are associated with Leber congenital amaurosis and juvenile retinitis pigmentosa. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	SPATA7
Immunogen:	Fusion protein of human SPATA7
Full name:	spermatogenesis associated 7
Synonyms :	HSD3; LCA3; HSD-3.1; HEL-S-296

SwissProt:	Q9P0W8
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer and human cervical cancer
IHC Recommend dilution:	50-200

