

兔抗 SPATA16 多克隆抗体

中文名称：兔抗 SPATA16 多克隆抗体

英文名称： Anti-SPATA16 rabbit polyclonal antibody

别名： SPGF6; NYD-SP12

抗原： SPATA16

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a testis-specific protein that belongs to the tetratricopeptide repeat-like superfamily. The encoded protein localizes to the Golgi apparatus and may play a role in spermatogenesis.
Applications:	ELISA, WB, IHC
Name of antibody:	SPATA16
Immunogen:	Fusion protein of human SPATA16
Full name:	spermatogenesis associated 16
Synonyms:	SPGF6; NYD-SP12
SwissProt:	Q9BXB7
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human esophagus cancer
IHC Recommend dilution:	150-300
WB Predicted band size:	65 kDa
WB Positive control:	Mouse brain tissue lysate
WB Recommended dilution:	500-2000



