

兔抗 SPATA18 多克隆抗体

中文名称：兔抗 SPATA18 多克隆抗体

英文名称： Anti-SPATA18 rabbit polyclonal antibody

别名： Mieap; SPETEX1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SPATA18

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	SPATA18 (spermatogenesis associated 18), also known as Mieap or SPETEX1, is a 538 amino acid protein that is thought to play a role in cell differentiation during spermatogenesis, particularly during development from late elongate spermatids to mature spermatozoa. Localizing to cytoplasm, SPATA18 is encoded by a gene that maps to human chromosome 4q12.
Applications:	ELISA, WB
Name of antibody:	SPATA18
Immunogen:	Synthetic peptide of human SPATA18
Full name:	spermatogenesis associated 18

Synonyms :	Mieap; SPETEX1
SwissProt:	Q8TC71
ELISA Recommended dilution:	2000-5000
WB Predicted band size:	61 kDa
WB Positive control:	Human normal stomach tissue
WB Recommended dilution:	500-2000

