

兔抗 LIMK2(Ab-283) 多克隆抗体

中文名称: 兔抗 LIMK2(Ab-283) 多克隆抗体

英文名称: Anti-LIMK2(Ab-283) rabbit polyclonal antibody

储存: 冷冻 (-20℃) 避光

抗原: LIMK2(Ab-283)

宿主: Rabbit

反应种属: Human Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: Unconjugate

技术规格

Background:	There are approximately 40 known eukaryotic LIM proteins, so named for the LIM domains they contain. LIM domains are highly conserved cysteine-rich structures containing 2 zinc fingers. Although zinc fingers usually function by binding to DNA or RNA, the LIM motif probably mediates protein-protein interactions. LIM kinase-1 and LIM kinase-2 belong to a small subfamily with a unique combination of 2 N-terminal LIM motifs and a C-terminal protein kinase domain.
Applications:	WB
Name of antibody:	LIMK2(Ab-283)
Immunogen:	Synthesized non-phosphopeptide derived from human LIMK2 around the phosphorylation site of serine 283 (R-R-S(p)-L-R).
Full name:	LIM domain kinase 2

Synonyms :	0
SwissProt:	P53671
WB Predicted band size:	72 kDa
WB Positive control:	NIH/3T3 cells treated with PMA (125ng/ml, 30mins) lysate
WB Recommended dilution:	500-3000

