

兔抗 SPDL1 多克隆抗体

中文名称：兔抗 SPDL1 多克隆抗体

英文名称： Anti-SPDL1 rabbit polyclonal antibody

别名： CCDC99

抗原： SPDL1

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Required for the localization of dynein and dynactin to the mitotic kinetochore. Dynein is believed to control the initial lateral interaction between the kinetochore and spindle microtubules and to facilitate the subsequent formation of end-on kinetochore-microtubule attachments mediated by the NDC80 complex. Also required for correct spindle orientation. Does not appear to be required for the removal of spindle assembly checkpoint (SAC) proteins from the kinetochore upon bipolar spindle attachment.
Applications:	ELISA, IHC
Name of antibody:	SPDL1

Immunogen:	Synthetic peptide of humanSPDL1
Full name:	spindle apparatus coiled-coil protein 1
Synonyms :	CCDC99
SwissProt:	Q96EA4
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer and Human esophagus cancer
IHC Recommend dilution:	25-100

