

## 兔抗 LIN28A 多克隆抗体

中文名称：兔抗 LIN28A 多克隆抗体

英文名称：Anti-LIN28A rabbit polyclonal antibody

别名：CSDD1; LIN28; LIN-28; ZCCHC1; lin-28A

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：LIN28A

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

**Background:**

This gene encodes a LIN-28 family RNA-binding protein that acts as a posttranscriptional regulator of genes involved in developmental timing and self-renewal in embryonic stem cells. The encoded protein functions through direct interaction with target mRNAs and by disrupting the maturation of certain miRNAs involved in embryonic development. This protein prevents the terminal processing of the LET7 family of microRNAs which are major regulators of cellular growth and differentiation. Aberrant expression of this gene is associated with cancer progression in multiple tissues.

<b>Applications:</b>	WB
<b>Name of antibody:</b>	LIN28A
<b>Immunogen:</b>	Synthetic peptide of human LIN28A
<b>Full name:</b>	lin-28 homolog A
<b>Synonyms:</b>	CSDD1; LIN28; LIN-28; ZCCHC1; lin-28A
<b>SwissProt:</b>	Q9H9Z2
<b>WB Predicted band size:</b>	23 kDa
<b>WB Positive control:</b>	HepG2 cells
<b>WB Recommended dilution:</b>	500-1000

