

兔抗 LIPE 多克隆抗体

中文名称：兔抗 LIPE 多克隆抗体

英文名称：Anti-LIPE rabbit polyclonal antibody

别名：HSL, LHS

抗原：LIPE

储存：冷冻（-20℃）

宿主：Rabbit

相关类别：一抗

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene has a long and a short form, generated by use of alternative translational start codons. The long form is expressed in steroidogenic tissues such as testis, where it converts cholesteryl esters to free cholesterol for steroid hormone production. The short form is expressed in adipose tissue, among others, where it hydrolyzes stored triglycerides to free fatty acids.
Applications:	ELISA, IHC
Name of antibody:	LIPE
Immunogen:	Fusion protein of human LIPE
Full name:	lipase, hormone-sensitive

Synonyms :	HSL, LHS
SwissProt:	Q05469
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human esophagus cancer
IHC Recommend dilution:	10-50

