

兔抗 SPINK6 多克隆抗体

中文名称：兔抗 SPINK6 多克隆抗体

英文名称： Anti-SPINK6 rabbit polyclonal antibody

别名： BUSI2; UNQ844

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： SPINK6

反应种属： Human

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	The protein encoded by this gene is a Kazal-type serine protease inhibitor that acts on kallikrein-related peptidases in the skin. Two transcript variants the same protein have been found for this gene.
Applications:	IHC
Name of antibody:	SPINK6
Immunogen:	Synthesized peptide derived from internal of human SPINK6.
Full name:	serine peptidase inhibitor, Kazal type 6
Synonyms :	BUSI2; UNQ844
SwissProt:	Q6UWN8

IHC positive control:	Human lung carcinoma tissue
IHC Recommend dilution:	50-100

