

## 兔抗 SPN 多克隆抗体

中文名称：兔抗 SPN 多克隆抗体

英文名称： Anti-SPN rabbit polyclonal antibody

别名： LSN; CD43; GALGP; GPL115

相关类别： 一抗

抗原： SPN

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a major sialoglycoprotein found on the surface of thymocytes, T lymphocytes, monocytes, granulocytes, and some B lymphocytes. It may be part of a physiologic ligand-receptor complex involved in T-cell activation. During T-cell activation, this protein is actively removed from the T-cell-APC (antigen-presenting cell) contact site, suggesting a negative regulatory role in adaptive immune response.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	SPN
<b>Immunogen:</b>	Fusion protein of human SPN

<b>Full name:</b>	sialophorin
<b>Synonyms :</b>	LSN; CD43; GALGP; GPL115
<b>SwissProt:</b>	P16150
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human tonsil and Human ovarian cancer
<b>IHC Recommend dilution:</b>	50-200

