

兔抗 SPOCK3 多克隆抗体

中文名称：兔抗 SPOCK3 多克隆抗体

英文名称： Anti-SPOCK3 rabbit polyclonal antibody

别名： SPARC/osteonectin, cwcv and kazal like domains proteoglycan 3; TES-3; TICN3; HSAJ1454

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SPOCK3

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of a novel family of calcium-binding proteoglycan proteins that contain thyroglobulin type-1 and Kazal-like domains. The encoded protein and may play a role in adult T-cell leukemia by inhibiting the activity of membrane-type matrix metalloproteinases. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.
Applications:	ELISA, WB, IHC
Name of antibody:	SPOCK3

Immunogen:	Fusion protein of human SPOCK3
Full name:	SPARC/osteonectin, cwcv and kazal like domains prot eoglycan 3
Synonyms:	TES-3; TICN3; HSAJ1454
SwissProt:	Q9BQ16
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	25-100
WB Predicted band size:	49 kDa
WB Positive control:	Human left kidney tissue,Human cerebrum tissue,Mo use brain tissue,HL-60 cell lysates
WB Recommended dilution:	200-1000



