

兔抗 LMO3 多克隆抗体

中文名称: 兔抗 LMO3 多克隆抗体

英文名称: Anti-LMO3 rabbit polyclonal antibody

别名: LIM domain only 3; RBTN3; RHOM3; RBTNL2; Rhom-3

抗原: LMO3

储存: 冷冻 (-20℃)

宿主: Rabbit

相关类别: 一抗

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	The protein encoded by this gene belongs to the rhombotin family of cysteine-rich LIM domain oncogenes. This gene is predominantly expressed in the brain. Related family members, LMO1 and LMO2 on chromosome 11, have been reported to be involved in chromosomal translocations in T-cell leukemia. Many alternatively spliced transcript variants have been found for this gene.
Applications:	ELISA, WB, IHC
Name of antibody:	LMO3
Immunogen:	Fusion protein of human LMO3
Full name:	LIM domain only 3

Synonyms:	RBTN3; RHOM3; RBTNL2; Rhom-3
SwissProt:	Q8TAP4
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	50-300
WB Predicted band size:	17 kDa
WB Positive control:	Human cerebella tissue, Rat heart tissue, Human heart tissue lysates
WB Recommended dilution:	200-1000



