

兔抗 SPON1 多克隆抗体

中文名称：兔抗 SPON1 多克隆抗体

英文名称： Anti-SPON1 rabbit polyclonal antibody

别名： f-spondin; VSGP/F-spondin

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SPON1

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	F-Spondin, also designated Spondin-1 or vascular smooth muscle growth-promoting factor, is a member of the subgroup of the Thrombospondin type 1 class molecules. F-Spondin is a secreted, extracellular matrix-attached protein which patterns axonal trajectories by promoting adhesion and outgrowth of commissural axons, in addition to inhibiting outgrowth of motor axons.
Applications:	ELISA, WB, IHC
Name of antibody:	SPON1
Immunogen:	Synthetic peptide of human SPON1
Full name:	spondin 1, extracellular matrix protein

Synonyms :	f-spondin; VSGP/F-spondin
SwissProt:	Q9HCB6
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	91 kDa
WB Positive control:	Mouse brain and lung tissue
WB Recommended dilution:	200-1000



