

小鼠抗 ANXA4 单克隆抗体

- 中文名称： 小鼠抗 ANXA4 单克隆抗体
- 英文名称： Anti-ANXA4 mouse monoclonal antibody
- 相关类别： 一抗
- 储 存： 冷冻（-20℃）
- 宿 主： Mouse
- 抗 原： ANXA4
- 反应种属： Human
- 标 记 物： Unconjugate
- 克隆类型： mouse monoclonal

技术规格

Background:

Annexin IV (ANX4) belongs to the annexin family of calcium-dependent phospholipid binding proteins. Although their functions are still not clearly defined, several members of the annexin family have been implicated in membrane-related events along exocytotic and endocytotic pathways. ANX4 has 45 to 59% identity with other members of its family and shares a similar size and exon-intron organization. Isolated from human placenta, ANX4 encodes a protein that has possible interactions with ATP, and has in vitro anticoagulant activity and also inhibits phospholipase A2 activity. ANX4 is almost exclusively expressed in epithelial cells. Several transcript variants encoding different isoforms have been found for this gene.

Applications:	ELISA, WB
Name of antibody:	ANXA4
Immunogen:	Fusion protein of human ANXA4
Full name:	annexin A4
Synonyms:	ANX4; P32.5; PIG28; PP4-X; ZAP36; PAP-II; HEL-S-274
SwissProt:	P09525
WB Predicted band size:	36 kDa
WB Positive control:	H520 whole cell lysate
WB Recommended dilution:	200-1000

