

## 兔抗 LMO4 多克隆抗体

中文名称：兔抗 LMO4 多克隆抗体

英文名称：Anti-LMO4 rabbit polyclonal antibody

相关类别：一抗

储 存：冷冻（-20℃）

宿 主：Rabbit

抗 原：LMO4

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a cysteine-rich protein that contains two LIM domains but lacks a DNA-binding homeodomain. The encoded protein may play a role as a transcriptional regulator or as an oncogene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	LMO4
<b>Immunogen:</b>	Fusion protein of human LMO4
<b>Full name:</b>	LIM domain only 4
<b>SwissProt:</b>	P61968
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human esophagus cancer
<b>IHC Recommend dilution:</b>	20-100

