

## 小鼠抗 ANXA10 单克隆抗体

中文名称： 小鼠抗 ANXA10 单克隆抗体

英文名称： Anti-ANXA10 mouse monoclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃）

宿 主： Mouse

抗 原： ANXA10

反应种属： Human

标 记 物： Unconjugate

克隆类型： mouse monoclonal、

### 技术规格

<b>Background:</b>	This gene encodes a member of the annexin family. Members of this calcium-dependent phospholipid-binding protein family play a role in the regulation of cellular growth and in signal transduction pathways. The function of this gene has not yet been determined.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	ANXA10
<b>Immunogen:</b>	Fusion protein of human ANXA10
<b>Full name:</b>	annexin A10
<b>Synonyms:</b>	ANX14
<b>SwissProt:</b>	Q9UJ72
<b>WB Predicted band size:</b>	37 kDa
<b>WB Positive control:</b>	HeLa whole cell lysate

WB Recommended dilution: 500-2000

