

### 兔抗 SPRY4 多克隆抗体

中文名称：兔抗 SPRY4 多克隆抗体

英文名称： Anti-SPRY4 rabbit polyclonal antibody

别名： HH17

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SPRY4

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of a family of cysteine- and proline-rich proteins. The encoded protein is an inhibitor of the receptor-transduced mitogen-activated protein kinase (MAPK) signaling pathway. Activity of this protein impairs the formation of active GTP-RAS. Nucleotide variation in this gene has been associated with hypogonadotropic hypogonadism 17 with or without anosmia. Alternative splicing results in a multiple transcript variants.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	SPRY4
<b>Immunogen:</b>	Full length fusion protein

<b>Full name:</b>	sprouty RTK signaling antagonist 4
<b>Synonyms :</b>	HH17
<b>SwissProt:</b>	Q9C004
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human lung cancer
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	33 kDa
<b>WB Positive control:</b>	K562 , A375 and 231 cell
<b>WB Recommended dilution:</b>	500-2000



