

小鼠抗 APOA4 单克隆抗体

中文名称： 小鼠抗 APOA4 单克隆抗体

英文名称： Anti-APOA4 mouse monoclonal antibody

别名： apolipoprotein A4

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： APOA4

反应种属： Human

标记物： Unconjugate

克隆类型： mouse monoclonal

技术规格

Background:	Apolipoprotein (apo) A-IV gene contains 3 exons separated by two introns. A sequence polymorphism has been identified in the 3'UTR of the third exon. The primary translation product is a 396-residue preprotein which after proteolytic processing is secreted its primary site of synthesis, the intestine, in association with chylomicron particles. Although its precise function is not known, apo A-IV is a potent activator of lecithin-cholesterol acyltransferase in vitro.
Applications:	ELISA, WB
Name of antibody:	APOA4
Immunogen:	Fusion protein of human APOA4

Full name:	apolipoprotein A4
SwissProt:	P06727
WB Predicted band size:	45 kDa
WB Positive control:	MCF-7 whole cell lysate
WB Recommended dilution:	200-1000

