

## 兔抗 SPSB2 多克隆抗体

中文名称：兔抗 SPSB2 多克隆抗体

英文名称： Anti-SPSB2 rabbit polyclonal antibody

别名： splA/ryanodine receptor domain and SOCS box containing 2; SSB2; GRCC9

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SPSB2

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of a subfamily of proteins containing a central SPRY (repeats in splA and RyR) domain and a C-terminal suppressor of cytokine signaling (SOCS) box. This protein plays a role in cell signaling. This gene is present in a gene-rich cluster on chromosome 12p13 in the vicinity of the CD4 antigen and triosephosphate isomerase genes. Alternative splicing results in multiple transcript variants.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	SPSB2
<b>Immunogen:</b>	Fusion protein of human SPSB2

<b>Full name:</b>	splA/ryanodine receptor domain and SOCS box containing 2
<b>Synonyms:</b>	SSB2; GRCC9
<b>SwissProt:</b>	Q99619
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer
<b>IHC Recommend dilution:</b>	50-100

