

兔抗 LPAR6 多克隆抗体

中文名称: 兔抗 LPAR6 多克隆抗体

英文名称: Anti-LPAR6 rabbit polyclonal antibody

别名: LAH3; P2Y5; ARWH1; HYPT8; P2RY5

相关类别: 一抗

储存: 冷冻 (-20℃)

抗原: LPAR6

宿主: Rabbit

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	The protein encoded by this gene belongs to the family of G-protein coupled receptors, that are preferentially activated by adenosine and uridine nucleotides. This gene aligns with an internal intron of the retinoblastoma susceptibility gene in the reverse orientation. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	LPAR6
Immunogen:	Synthetic peptide of human LPAR6
Full name:	lysophosphatidic acid receptor 6
Synonyms:	LAH3; P2Y5; ARWH1; HYPT8; P2RY5

SwissProt:	P43657
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer and human ovarian cancer
IHC Recommend dilution:	10-50



