

小鼠抗 ARHGAP5 单克隆抗体

中文名称： 小鼠抗 ARHGAP5 单克隆抗体

英文名称： Anti-ARHGAP5 mouse monoclonal antibody

别名： Rho GTPase activating protein 5; GFI2; p190-B; RhoGAP5; p190BRhoGAP

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： ARHGAP5

反应种属： Human

标记物： Unconjugate

克隆类型： mouse monoclonal

技术规格

Background:	Rho GTPase activating protein 5 negatively regulates RHO GTPases, a family which may mediate cytoskeleton changes by stimulating the hydrolysis of bound GTP. Two transcript variants encoding different isoforms have been found for this gene.
Applications:	WB
Name of antibody:	ARHGAP5
Immunogen:	Fusion protein of human ARHGAP5
Full name:	Rho GTPase activating protein 5
Synonyms:	GFI2; p190-B; RhoGAP5; p190BRhoGAP
SwissProt:	Q13017



WB Predicted band size:	172 KD
WB Positive control:	Hela, NIH/3T3, Jurkat, K562, MCF7, A431 cell lysates
WB Recommended dilution:	200-1000