

兔抗 SRC(Phospho-Ser75) 多克隆抗体

中文名称：兔抗 SRC(Phospho-Ser75) 多克隆抗体

英文名称： Anti-SRC(Phospho-Ser75) rabbit polyclonal antibody

别名： ASV; SRC1; c-SRC; p60-Src

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： SRC(Phospho-Ser75)

反应种属： Human Mouse Rat

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	This gene is highly similar to the v-src gene of Rous sarcoma virus. This proto-oncogene may play a role in the regulation of embryonic development and cell growth. The protein encoded by this gene is a tyrosine-protein kinase whose activity can be inhibited by phosphorylation by c-SRC kinase. Mutations in this gene could be involved in the malignant progression of colon cancer. Two transcript variants encoding the same protein have been found for this gene.
Applications:	WB

Name of antibody:	SRC(Phospho-Ser75)
Immunogen:	Peptide sequence around phosphorylation site of serine 75 (V-T-S(p)-P-Q) derived from Human Src.
Full name:	SRC proto-oncogene, non-receptor tyrosine kinase
Synonyms :	ASV; SRC1; c-SRC; p60-Src
SwissProt:	P12931
WB Predicted band size:	60 kDa
WB Positive control:	COLO205lysates
WB Recommended dilution:	500-1000

