

小鼠抗 ARHGEF18 单克隆抗体

中文名称： 小鼠抗 ARHGEF18 单克隆抗体

英文名称： Anti-ARHGEF18 mouse monoclonal antibody

别名： Rho/Rac guanine nucleotide exchange factor 18; RP78; P114-RhoGEF

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： ARHGEF18

反应种属： Human

标记物： Unconjugate

克隆类型： mouse monoclonal

技术规格

Background:

Rho GTPases are GTP binding proteins that regulate a wide spectrum of cellular functions. These cellular processes include cytoskeletal rearrangements, gene transcription, cell growth and motility. Activation of Rho GTPases is under the direct control of guanine nucleotide exchange factors (GEFs). The protein encoded by this gene is a guanine nucleotide exchange factor and belongs to the Rho GTPase GFE family. Family members share a common feature, a Dbl (DH) homology domain followed by a pleckstrin (PH) homology domain. Alternatively spliced transcript variants encoding different isoforms have been identified.

Applications:	WB
Name of antibody:	ARHGEF18
Immunogen:	Fusion protein of human ARHGEF18
Full name:	Rho/Rac guanine nucleotide exchange factor 18
Synonyms:	RP78; P114-RhoGEF
SwissProt:	Q6ZSZ5
WB Predicted band size:	131 KD
WB Positive control:	293T, A431, Hela, HepG2 cell lysates
WB Recommended dilution:	500-2000