

## 兔抗 SRC(Phospho-Tyr418) 多克隆抗体

中文名称：兔抗 SRC(Phospho-Tyr418) 多克隆抗体

英文名称：Anti-SRC(Phospho-Tyr418) rabbit polyclonal antibody

别名：ASV; SRC1; c-SRC; p60-Src

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：SRC(Phospho-Tyr418)

反应种属：Human Mouse Rat

标记物：Unconjugate

克隆类型：Unconjugate

### 技术规格

<b>Background:</b>	This gene is highly similar to the v-src gene of Rous sarcoma virus. This proto-oncogene may play a role in the regulation of embryonic development and cell growth. The protein encoded by this gene is a tyrosine-protein kinase whose activity can be inhibited by phosphorylation by c-SRC kinase. Mutations in this gene could be involved in the malignant progression of colon cancer. Two transcript variants encoding the same protein have been found for this gene.
<b>Applications:</b>	WB, IHC

<b>Name of antibody:</b>	SRC(Phospho-Tyr418)
<b>Immunogen:</b>	Peptide sequence around phosphorylation site of tyrosine 418(N-E-Y(p)-T-A)derived from Human Src.
<b>Full name:</b>	SRC proto-oncogene, non-receptor tyrosine kinase
<b>Synonyms :</b>	ASV; SRC1; c-SRC; p60-Src
<b>SwissProt:</b>	P12931
<b>IHC positive control:</b>	Human breast carcinoma tissue
<b>IHC Recommend dilution:</b>	50-100
<b>WB Predicted band size:</b>	60 kDa
<b>WB Positive control:</b>	COLO205 cells lysates
<b>WB Recommended dilution:</b>	500-1000



