

兔抗 SRCIN1 多克隆抗体

中文名称： 兔抗 SRCIN1 多克隆抗体

英文名称： Anti-SRCIN1 rabbit polyclonal antibody

别名： P140; SNIP

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： SRCIN1

反应种属： Human

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	Acts as a negative regulator of SRC by activating CSK which inhibits SRC activity and downstream signaling, leading to impaired cell spreading and migration. Regulates dendritic spine morphology. Involved in calcium-dependent exocytosis. May play a role in neurotransmitter release or synapse maintenance.
Applications:	WB
Name of antibody:	SRCIN1
Immunogen:	Synthesized peptide derived from internal of human SNIP.
Full name:	SRC kinase signaling inhibitor 1

Synonyms :	P140; SNIP
SwissProt:	Q9C0H9
WB Predicted band size:	112 kDa
WB Positive control:	HepG2 and Hela cells lysates
WB Recommended dilution:	500-3000

