

小鼠抗 ARPC5 单克隆抗体

中文名称: 小鼠抗 ARPC5 单克隆抗体

英文名称: Anti-ARPC5 mouse monoclonal antibody

别名: ARC16; dJ127C7.3; p16-Arc

抗原: ARPC5

储存: 冷冻 (-20℃) 避光

宿主: Mouse

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: mouse monoclonal

技术规格

Background:	This gene encodes one of seven subunits of the human Arp2/3 protein complex. The Arp2/3 protein complex has been implicated in the control of actin polymerization in cells and has been conserved through evolution. The exact role of the protein encoded by this gene, the p16 subunit, has yet to be determined. Alternatively spliced transcript variants encoding different isoforms have been observed for this gene. [provided by RefSeq, Jul 2013]
Applications:	WB, IHC
Name of antibody:	ARPC5
Immunogen:	Fusion protein of human ARPC5

Full name:	actin related protein 2/3 complex, subunit 5, 16kDa (ARPC5), transcript variant 1
Synonyms:	ARC16; dJ127C7.3; p16-Arc
SwissProt:	O15511
IHC positive control:	adenocarcinoma of human ovary tissue and adenocarcinoma of human colon tissue
IHC Recommend dilution:	30-150
WB Predicted band size:	16 kDa
WB Positive control:	A549, Jurkat, MDCK, MCF-7 cell lysates
WB Recommended dilution:	200-1000