

兔抗 LRATD2 多克隆抗体

- 中文名称：兔抗 LRATD2 多克隆抗体
- 英文名称：Anti-LRATD2 rabbit polyclonal antibody
- 别名：LRAT domain containing 2; NSE2; FAM84B; BCMP101
- 相关类别：一抗
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 抗原：LRATD2
- 反应种属：Human
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

技术规格

Background:	FAM84B (family with sequence similarity 84, member B), also known as NSE2 or BCMP101, is a 310 amino acid protein that is expressed in esophageal squamous cell carcinomas, suggesting a role in tumor development and metastasis. The gene encoding FAM84B maps to human chromosome 8, which consists of nearly 146 million base pairs, houses more than 800 genes and is associated with a variety of diseases and malignancies. Schizophrenia, bipolar disorder, trisomy 8, Pfeiffer syndrome, congenital hypothyroidism, Waardenburg syndrome and some leukemias and lymphomas are thought to occur as a result of defects in specific genes that map to chr
--------------------	--

	omosome 8.
Applications:	ELISA, WB, IHC
Name of antibody:	LRATD2
Immunogen:	Synthetic peptide of human LRATD2
Full name:	LRAT domain containing 2
Synonyms:	NSE2; FAM84B; BCMP101
SwissProt:	Q96KN1
IHC positive control:	Human esophagus cancer and Human thyroid cancer
IHC Recommend dilution:	50-300
WB Predicted band size:	34 kDa
WB Positive control:	HUVEC and PC3 cell lysates
WB Recommended dilution:	500-2000



