

兔抗 LRIF1 多克隆抗体

中文名称：兔抗 LRIF1 多克隆抗体

英文名称：Anti-LRIF1 rabbit polyclonal antibody

别名：ligand dependent nuclear receptor interacting factor 1; RIF1; C1orf103

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：LRIF1

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

| | |
|--------------------------|---|
| Background: | Represses the ligand-induced transcriptional activity of retinoic acid receptor alpha (RARA). This repression may occur through direct recruitment of histone deacetylases. |
| Applications: | ELISA, IHC |
| Name of antibody: | LRIF1 |
| Immunogen: | Fusion protein of human LRIF1 |
| Full name: | ligand dependent nuclear receptor interacting factor 1 |
| Synonyms: | RIF1; C1orf103 |
| SwissProt: | Q5T3J3 |

| | |
|------------------------------------|-----------------------|
| ELISA Recommended dilution: | 5000-10000 |
| IHC positive control: | Human cervical cancer |
| IHC Recommend dilution: | 50-200 |

