

## 兔抗 SRGN 多克隆抗体

中文名称：兔抗 SRGN 多克隆抗体

英文名称： Anti-SRGN rabbit polyclonal antibody

别名： PPG; PRG; PRG1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SRGN

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a protein best known as a hematopoietic cell granule proteoglycan. Proteoglycans stored in the secretory granules of many hematopoietic cells also contain a protease-resistant peptide core, which may be important for neutralizing hydrolytic enzymes. This encoded protein was found to be associated with the macromolecular complex of granzymes and perforin, which may serve as a mediator of granule-mediated apoptosis.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	SRGN

<b>Immunogen:</b>	Fusion protein of human SRGN
<b>Full name:</b>	serglycin
<b>Synonyms :</b>	PPG; PRG; PRG1
<b>SwissProt:</b>	P10124
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human liver cancer
<b>IHC Recommend dilution:</b>	25-100

