

小鼠抗 AZGP1 单克隆抗体

中文名称：小鼠抗 AZGP1 单克隆抗体

英文名称：Anti-AZGP1 mouse monoclonal antibody

别名：alpha-2-glycoprotein 1, zinc-binding; ZAG; ZA2G

相关类别：一抗

储存：冷冻（-20℃）

宿主：Mouse

抗原：AZGP1

反应种属：Human

标记物：Unconjugate

克隆类型：mouse monoclonal

技术规格

Background:	Stimulates lipid degradation in adipocytes and causes the extensive fat losses associated with some advanced cancers. May bind polyunsaturated fatty acids.
Applications:	WB, IHC
Name of antibody:	AZGP1
Immunogen:	Fusion protein of human AZGP1
Full name:	alpha-2-glycoprotein 1, zinc-binding
Synonyms:	ZAG; ZA2G
SwissProt:	P25311
IHC positive control:	Human liver cancer and Human liver tissue; Human breast tissue
IHC Recommend dilution:	500-2000

WB Predicted band size:	34 KD
WB Positive control:	Human liver tissue, Human heart tissue, Human lung tissue lysates
WB Recommended dilution:	500-2000