

## 小鼠抗 BECN1 单克隆抗体

中文名称： 小鼠抗 BECN1 单克隆抗体

英文名称： Anti-BECN1 mouse monoclonal antibody

别名： beclin 1; ATG6; VPS30; beclin1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： BECN1

反应种属： Human, Mouse, Rat, Monkey, Hamster

标记物： Unconjugate

克隆类型： Mouse Monoclonal

### 技术规格

<b>Background:</b>	This gene encodes a protein that regulates autophagy, a catabolic process of degradation induced by starvation. The encoded protein is a component of the phosphatidylinositol-3-kinase (PI3K) complex which mediates vesicle trafficking processes. This protein is thought to play a role in multiple cellular processes, including tumorigenesis, neurodegeneration and apoptosis. Alternative splicing results in multiple transcript variants.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	BECN1
<b>Immunogen:</b>	Fusion protein of human BECN1

<b>Full name:</b>	beclin 1
<b>Synonyms:</b>	ATG6; VPS30; beclin1
<b>SwissProt:</b>	Q14457
<b>WB Predicted band size:</b>	52 kDa
<b>WB Positive control:</b>	Jurkat cell, Rat brain tissue, CHO-K1, COS7, C6, 3T3 and Hela cell lysates
<b>WB Recommended dilution:</b>	1000-5000

