

兔抗 SSBP1 多克隆抗体

中文名称：兔抗 SSBP1 多克隆抗体

英文名称： Anti-SSBP1 rabbit polyclonal antibody

别名： SSBP; mtSSB; Mt-SSB; SOSS-B1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SSBP1

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	SSBP1 is a housekeeping gene involved in mitochondrial biogenesis (Tiranti et al., 1995 [PubMed 7789991]). It is also a subunit of a single-stranded DNA (ssDNA)-binding complex involved in the maintenance of genome stability.
Applications:	ELISA, IHC
Name of antibody:	SSBP1
Immunogen:	Fusion protein of human SSBP1
Full name:	single-stranded DNA binding protein 1, mitochondrial
Synonyms :	SSBP; mtSSB; Mt-SSB; SOSS-B1

SwissProt:	Q04837
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer and human esophagus cancer
IHC Recommend dilution:	25-100

