

## 兔抗 SSTR2 多克隆抗体

中文名称：兔抗 SSTR2 多克隆抗体

英文名称：Anti-SSTR2 rabbit polyclonal antibody

相关类别：一抗

储 存：冷冻（-20℃）

宿 主：Rabbit

抗 原：SSTR2

反应种属：Human, Mouse, Rat

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Somatostatin acts at many sites to inhibit the release of many hormones and other secretory proteins. The biologic effects of somatostatin are probably mediated by a family of G protein-coupled receptors that are expressed in a tissue-specific manner. SSTR2 is a member of the superfamily of receptors having seven transmembrane segments and is expressed in highest levels in cerebrum and kidney.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	SSTR2
<b>Immunogen:</b>	Synthetic peptide of human SSTR2
<b>Full name:</b>	somatostatin receptor 2
<b>SwissProt:</b>	P30874

<b>ELISA Recommended dilution:</b>	1000-5000
<b>WB Predicted band size:</b>	41 kDa
<b>WB Positive control:</b>	Mouse thymus tissue
<b>WB Recommended dilution:</b>	200-1000

