

兔抗 SST 多克隆抗体

中文名称：兔抗 SST 多克隆抗体

英文名称： Anti-SST rabbit polyclonal antibody

别名： somatostatin; SMST

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SST

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

The hormone somatostatin has active 14 aa and 28 aa forms that are produced by alternate cleavage of the single preproprotein encoded by this gene. Somatostatin is expressed throughout the body and inhibits the release of numerous secondary hormones by binding to high-affinity G-protein-coupled somatostatin receptors. This hormone is an important regulator of the endocrine system through its interactions with pituitary growth hormone, thyroid stimulating hormone, and most hormones of the gastrointestinal tract. Somatostatin also affects rates of neurotransmission in the central nervous

	system and proliferation of both normal and tumorigenic cells.
Applications:	ELISA, WB, IHC
Name of antibody:	SST
Immunogen:	Synthetic peptide of human SST
Full name:	somatostatin
Synonyms:	SMST
SwissProt:	P61278
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human brain
IHC Recommend dilution:	25-50
WB Predicted band size:	13 kDa
WB Positive control:	A549 and Hela cell lysates
WB Recommended dilution:	200-1000



