

小鼠抗 BTLA 单克隆抗体

中文名称: 小鼠抗 BTLA 单克隆抗体

英文名称: Anti-BTLA mouse monoclonal antibody

别名: BTLA1; CD272

储存: 冷冻 (-20℃) 避光

抗原: BTLA

宿主: Mouse

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: mouse monoclonal

技术规格

Background:	BTLA expression is induced during activation of T cells, and BTLA remains expressed on Th1 cells but not Th2 cells. Like PD1 (MIM 600244) and CTLA4 (MIM 123890), BTLA interacts with a B7 homolog, B7H4.[supplied by OMIM].
Applications:	WB, IF, FC
Name of antibody:	BTLA
Immunogen:	Fusion protein of human BTLA
Full name:	B and T lymphocyte associated (BTLA), transcript variant 1
Synonyms:	BTLA1; CD272
SwissProt:	Q7Z6A9
WB Predicted band size:	33 kDa

WB Positive control:	Human testis tissue, Human uterus tissue, Human breast tissue, Human brain tissue, Human liver tissue, Human ovary tissue lysates
WB Recommended dilution:	500-2000