

## 兔抗 SSX5 多克隆抗体

中文名称：兔抗 SSX5 多克隆抗体

英文名称： Anti-SSX5 rabbit ployclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃）

宿 主： Rabbit

抗 原： SSX5

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The product of this gene belongs to the family of highly homologous synovial sarcoma X (SSX) breakpoint proteins. These proteins may function as transcriptional repressors. They are also capable of eliciting spontaneous humoral and cellular immune responses in cancer patients, and are potentially useful targets in cancer vaccine-based immunotherapy. While some of the related SSX genes are involved in t(X;18)(p11.2;q11.2) translocations that are characteristically found in all synovial sarcomas, this gene does not appear to be involved in such translocations.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	SSX5

<b>Immunogen:</b>	Synthetic peptide of human SSX5
<b>Full name:</b>	synovial sarcoma, X breakpoint 5
<b>SwissProt:</b>	O60225
<b>ELISA Recommended dilution:</b>	1000-2000
<b>WB Predicted band size:</b>	22 kDa
<b>WB Positive control:</b>	A431, Jurkat and Raji cells
<b>WB Recommended dilution:</b>	200-1000

