

## 小鼠抗 C2CD6 单克隆抗体

中文名称：小鼠抗 C2CD6 单克隆抗体

英文名称：Anti-C2CD6 mouse monoclonal antibody

别名：C2 calcium dependent domain containing 6; ALS2CR11

相关类别：一抗

储存：冷冻（-20℃）

宿主：Mouse

抗原：C2CD6

反应种属：Human

标记物：Unconjugate

克隆类型：mouse monoclonal

### 技术规格

<b>Background:</b>	An autosomal recessive form of juvenile amyotrophic lateral sclerosis was originally mapped to a region of chromosome 2 that includes this gene. The encoded protein contains a calcium-dependent membrane targeting C2 domain. This domain is often found in proteins that are involved in membrane trafficking and signal transduction.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	C2CD6
<b>Immunogen:</b>	Fusion protein of human C2CD6
<b>Full name:</b>	C2 calcium dependent domain containing 6
<b>Synonyms:</b>	ALS2CR11

<b>SwissProt:</b>	Q53TS8
<b>WB Predicted band size:</b>	71 KD
<b>WB Positive control:</b>	A375 cell lysate
<b>WB Recommended dilution:</b>	200-1000