

## 兔抗 STAB1 多克隆抗体

中文名称：兔抗 STAB1 多克隆抗体

英文名称：Anti-STAB1 rabbit polyclonal antibody

别名：FEX1; FEEL-1; FELE-1; STAB-1; CLEVER-1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：STAB

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a large, transmembrane receptor protein which may function in angiogenesis, lymphocyte homing, cell adhesion, or receptor scavenging. The protein contains 7 fasciclin, 16 epidermal growth factor (EGF)-like, and 2 laminin-type EGF-like domains as well as a C-type lectin-like hyaluronan-binding Link module. The protein is primarily expressed on sinusoidal endothelial cells of liver, spleen, and lymph node.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	STAB1

<b>Immunogen:</b>	Synthetic peptide of human STAB1
<b>Full name:</b>	stabilin 1
<b>Synonyms :</b>	FEX1; FEEL-1; FELE-1; STAB-1; CLEVER-1
<b>SwissProt:</b>	Q9NY15
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human liver cancer
<b>IHC Recommend dilution:</b>	25-100

