

## 兔抗 STAC2 多克隆抗体

- 中文名称: 兔抗 STAC2 多克隆抗体
- 英文名称: Anti-STAC2 rabbit polyclonal antibody
- 别名: 24b2; 24b2/STAC2
- 相关类别: 一抗
- 储存: 冷冻(-20℃)避光
- 宿 主: Rabbit
- 抗 原: STAC2
- 反应种属: Human
- 标记物: Unconjugate
- 克隆类型: Unconjugate

## 技术规格

P	
	STAC2 (SH3 and cysteine-rich domain-containin
Background:	g protein 2), also known as 24b2/STAC2, is a 4
	11 amino acid protein that contains one SH3 d
	omain and one phorbol-ester/DAG-type zinc fin
	ger. Encoding two alternatively spliced isoforms
	, the STAC2 gene contains over 15 thousand b
	ases and maps to human chromosome 17q12.
	Human chromosome 17 comprises over 2.5% o
	f the human genome and encodes over 1,200
	genes. Two key tumor suppressor genes are as
	sociated with chromosome 17, namely, p53 and



	BRCA1. Tumor suppressor p53 is necessary for maintenance of cellular genetic integrity by mo derating cell fate through DNA repair versus ce II death. Malfunction or loss of p53 expression is associated with malignant cell growth and Li -Fraumeni syndrome. Like p53, BRCA1 is directl y involved in DNA repair, though specifically it is recognized as a genetic determinant of early onset breast cancer and predisposition to canc ers of the ovary, colon, prostate gland and fall opian tubes.
Applications:	WB
Name of antibody:	STAC2
Immunogen:	Synthesized peptide derived from internal of h uman STAC2.
Full name:	SH3 and cysteine rich domain 2
Synonyms:	24b2; 24b2/STAC2
SwissProt:	Q6ZMT1
WB Predicted band size:	45 kDa
WB Positive control:	K562 cells lysate
WB Recommended dilution:	500-3000



