

兔抗 STAG1 多克隆抗体

中文名称：兔抗 STAG1 多克隆抗体

英文名称：Anti-STAG1 rabbit polyclonal antibody

别名：stromal antigen 1; SA1; MRD47; SCC3A

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：STAG1

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene is a member of the SCC3 family and is expressed in the nucleus. It encodes a component of cohesin, a multisubunit protein complex that provides sister chromatid cohesion along the length of a chromosome from DNA replication through prophase and prometaphase, after which it is dissociated in preparation for segregation during anaphase.
Applications:	ELISA, IHC
Name of antibody:	STAG1
Immunogen:	Fusion protein of human STAG1
Full name:	stromal antigen 1

Synonyms:	SA1; MRD47; SCC3A
SwissProt:	Q8WVM7
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer and Human tonsil
IHC Recommend dilution:	50-200



