

## 小鼠抗 CALU 单克隆抗体

中文名称： 小鼠抗 CALU 单克隆抗体

英文名称： Anti-CALU mouse monoclonal antibody

别名： calumenin

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： CALU

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： Mouse Monoclonal

### 技术规格

<b>Background:</b>	The product of this gene is a calcium-binding protein localized in the endoplasmic reticulum (ER) and it is involved in such ER functions as protein folding and sorting. This protein belongs to a family of multiple EF-hand proteins (CERC) that include reticulocalbin, ERC-55, and Cab45 and the product of this gene. Alternatively spliced transcript variants encoding different isoforms have been identified.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	CALU

<b>Immunogen:</b>	Fusion protein of human CALU
<b>Full name:</b>	calumenin
<b>SwissProt:</b>	O43852
<b>WB Predicted band size:</b>	37 kDa
<b>WB Positive control:</b>	SW480, mouse intestine and A431 cell lysates
<b>WB Recommended dilution:</b>	1000-5000

