

## 兔抗 STAG3 多克隆抗体

中文名称：兔抗 STAG3 多克隆抗体

英文名称： Anti-STAG3 rabbit polyclonal antibody

抗 原： STAG3

储 存： 冷冻（-20℃） 避光

宿 主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

### 技术规格

<b>Background:</b>	The protein encoded by this gene is expressed in the nucleus and is a subunit of the cohesin complex which regulates the cohesion of sister chromatids during cell division. A mutation in this gene is associated with premature ovarian failure. Alternate splicing results in multiple transcript variants encoding distinct isoforms. This gene has multiple pseudogenes.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	STAG3
<b>Immunogen:</b>	Synthesized peptide derived from C-terminal of human STAG3.
<b>Full name:</b>	stromal antigen 3
<b>Synonyms :</b>	0

<b>SwissProt:</b>	Q9UJ98
<b>WB Predicted band size:</b>	139 kDa
<b>WB Positive control:</b>	Jurkat cells lysate
<b>WB Recommended dilution:</b>	500-3000

