

## 小鼠抗 CAMK4 单克隆抗体

中文名称： 小鼠抗 CAMK4 单克隆抗体

英文名称： Anti-CAMK4 mouse monoclonal antibody

别名： calcium/calmodulin dependent protein kinase IV; caMK; CaMK IV; CaMK-GR

相关类别： 一抗

储存： 冷冻 (-20℃)

宿主： Mouse

抗原： CAMK4

反应种属： Human

标记物： Unconjugate

克隆类型： Mouse Monoclonal

### 技术规格

<b>Background:</b>	The product of this gene belongs to the serine/threonine protein kinase family, and to the Ca(2+)/calmodulin-dependent protein kinase subfamily. This enzyme is a multifunctional serine/threonine protein kinase with limited tissue distribution, that has been implicated in transcriptional regulation in lymphocytes, neurons and male germ cells.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	CAMK4
<b>Immunogen:</b>	Fusion protein of human CAMK4
<b>Full name:</b>	calcium/calmodulin dependent protein kinase IV
<b>Synonyms:</b>	caMK; CaMK IV; CaMK-GR

<b>SwissProt:</b>	Q16566
<b>WB Predicted band size:</b>	52 kDa
<b>WB Positive control:</b>	Romas and Jurkat cell lysates
<b>WB Recommended dilution:</b>	1000-5000

