

小鼠抗 CAPZA1 单克隆抗体

中文名称： 小鼠抗 CAPZA1 单克隆抗体

英文名称： Anti-CAPZA1 mouse monoclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃）

宿 主： Mouse

抗 原： CAPZA1

反应种属： Human

标 记 物： Unconjugate

克隆类型： mouse monoclonal

技术规格

Background:	CAPZA1 is a member of the F-actin capping protein alpha subunit family. This gene encodes the alpha subunit of the barbed-end actin binding protein. The protein regulates growth of the actin filament by capping the barbed end of growing actin filaments.
Applications:	ELISA, WB
Name of antibody:	CAPZA1
Immunogen:	Fusion protein of human CAPZA1
Full name:	capping actin protein of muscle Z-line alpha subunit 1
Synonyms:	CAPZ; CAZ1; CAPP1
SwissProt:	P52907
WB Predicted band size:	33 kDa
WB Positive control:	HT-29, HepG2, 293T, Jurkat and HeLa cell lysates

WB Recommended dilution: 500-2000

