

小鼠抗 CAPZA2 单克隆抗体

中文名称： 小鼠抗 CAPZA2 单克隆抗体

英文名称： Anti-CAPZA2 mouse monoclonal antibody

相关类别： 一抗

储 存： 冷冻 (-20℃)

宿 主： Mouse

抗 原： CAPZA2

反应种属： Human

标 记 物： Unconjugate

克隆类型： mouse monoclonal

技术规格

Background:	The protein encoded by this gene is a member of the F-actin capping protein alpha subunit family. It is the alpha subunit of the barbed-end actin binding protein Cap Z. By capping the barbed end of actin filaments, Cap Z regulates the growth of the actin filaments at the barbed end.
Applications:	ELISA, WB
Name of antibody:	CAPZA2
Immunogen:	Fusion protein of human CAPZA2
Full name:	capping actin protein of muscle Z-line alpha subunit 2
Synonyms:	CAPZ; CAPPA2; CapZ alpha-2
SwissProt:	P47755

WB Predicted band size:	33 kDa
WB Positive control:	SW480 whole cell lysate
WB Recommended dilution:	500-2000

KDa

66 —
47 —
42 —
33 —
26 —

