

兔抗 STARD10 多克隆抗体

中文名称：兔抗 STARD10 多克隆抗体

英文名称：Anti-STARD10 rabbit polyclonal antibody

别名：PCTP2; CGI-52; NY-CO-28; SDCCAG28

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：STARD10

反应种属：Human Mouse Rat

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	May play metabolic roles in sperm maturation or fertilization. By similarity. Phospholipid transfer protein that preferentially selects lipid species containing a palmitoyl or stearoyl chain on the sn-1 and an unsaturated fatty acyl chain (18:1 or 18:2) on the sn-2 position. Able to transfer phosphatidylcholine (PC) and phosphatidylethanolamine (PE) between membranes.
Applications:	WB
Name of antibody:	STARD10
Immunogen:	Synthesized peptide derived from N-terminal of human STARD10.

Full name:	StAR-related lipid transfer (START) domain containing 10
Synonyms :	PCTP2; CGI-52; NY-CO-28; SDCCAG28
SwissProt:	Q9Y365
WB Predicted band size:	33 kDa
WB Positive control:	COLO205 cells lysate
WB Recommended dilution:	500-3000

